

iPaxera

Mobile Viewer App



High performance, intuitive app to view DICOM studies anywhere, anytime from your tablet or smartphone

- Query / retrieve DICOM studies with any PACS system
- Display and navigate through images with interactive slider
- Control WW/WL with single finger drag
- Display Cine (multi-frame) images
- Seamless integration with PaxeraUtlima server
- View, write, dictate and send reports
- Image manipulations including zoom, pan, rotation, and distance measurement
- Connect other PACS via Q/R and WADO
- Import offline studies from iTunes™ and box™
- Navigate images instantly with unique ultra-fast streaming technology



Navigate studies



DICOM



Dictate Reports



Display Cine



Image Sharing



Image Manipulation

